Vakuumtechnik im System
Technology for Vacuum Systems
Technologie du vide

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Safety information!

NOTICE

Use the solenoid operated in-line valve for vacuum control by connecting it to a vacuum controller.

☞ Connect the solenoid operated valve only to a suitable vacuum controller (see "Technical data").

The **solenoid operated valve 20636668** is designed to be connected to the VACUUBRAND BUS system VACUU•BUS.

The **solenoid operated valve 20637659** is to be used only in combination with the vacuum controllers CVC 2II, CVC 2000 (only as venting valve), or VNC 1, or with the Vacuum Management System VMS Module A.

Remove all packing material from the packing box. Remove the product from its packing-box and retain all packaging for inspection until the equipment is inspected and tested. Inspect the equipment promptly and carefully. If the equipment is damaged, notify the supplier and the carrier in writing within three days. State the item number of the product together with the order number and the supplier's invoice number.

Use the equipment **for the intended use only**, i.e. as in-line valve in VACUU•LAN® Local Area Vacuum Networks or in vacuum systems. Do not use at overpressure. Do not use to distribute fluids or solids!

Comply with all applicable and relevant safety requirements (regulations and guidelines). **Implement the required actions and adopt suitable safety measures.**

Ensure that **maintenance** is done only by suitably trained and supervised technicians. As laid down in the statutory regulations (occupational, health and safety regulations and regulations for environmental protection), equipment returned to the manufacturer can be repaired only under certain conditions.

DANGER

➤ Adopt suitable measures to prevent the release of **dangerous, toxic, explosive, corrosive, noxious or polluting fluids, vapours and gases**.

➤ **Attention:** The equipment may be contaminated with chemicals, which have been processed during operation. Ensure that the equipment is completely decontaminated before maintenance commences. Take adequate precautions to protect people from the effects of dangerous substances if contamination has occurred. Use appropriate protective clothing, safety goggles and protective gloves.

WARNING

☞ Use only **OEM spare parts and accessories**. Otherwise safety and performance of the equipment might be reduced.

☞ **Do not use the equipment if it is damaged.**

☞ Check that the system is mechanically stable and that all fittings are secure.

☞ Do not disassemble any valve if under vacuum!

CAUTION

• Check that the power source and the equipment match in voltage, phase, and frequency (see "Technical data").

• Ensure that the materials of the wetted parts are compatible with the media, see "Technical data".

• **If idle, the solenoid operated valve is closed.**

Valve 2063668: The valve is open if energized and if Bus signal is applied.

Valve 20637659: If energized, the valve is open.

• Comply with max. permissible temperatures and pressures (see "Technical data").

• Check regularly that the connection between valve and vacuum system is leak-tight, especially prior to initial operation.

Assembly:

➤ Do not insert a sealing ring or a non-return valve.

➤ Tighten the **union nut** hand-tight; retighten after some days if necessary.

Technical data

Type	20636668	20637659
Maximum permissible range of supply voltage	24 VDC $\pm 10\%$	
Power draw (operation)	6 W	
Current draw	approx. 0.22 A	
Electrical connection	connector for VACUUBRAND BUS system VACUU•BUS	male connector 3-pole MAS 30 UM
Operating cycles per minute	maximum 50	
Valve currentless	valve closed	
Nominal diameter	DN 6	
Max. permissible operating pressure	1.5 bar	
Connection	union nut M35 x 1.5	
Max. permissible ambient temperature	40 °C	
Max. permissible gas temperature	80 °C	
Degree of protection IEC 529	IP 65	
Leak rate	$1 \cdot 10^{-2}$ mbar·l/s	
Dimensions L x W x H (without cable)	60 mm x 41 mm x 68 mm	
Cable length	2 m	
Weight (without cable)	0.21 kg	

Components	Wetted materials
Valve body	PVDF
Seal material	PTFE / fluorine elastomer

We reserve the right for technical modifications without prior notice!

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DECLARATION OF CONFORMITY – China RoHS 2

VACUUBRAND GMBH + CO KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(+VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Protection Use Period (EPUP)

EPUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The environmental Protection Use Period for VACUUBRAND products is 40 years.

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This table is created according to SJ/T 11363-2006.

MATERIAL CONTENT DECLARATION FOR VACUUBRAND PRODUCTS							
部件名称 Part name	有毒有害物质或元素 Hazardous substances						环保期限标识 EPUP
	铅 Pb	汞 Hg	镉 Cd	六价铬 Cr(+VI)	多溴联苯 PBB	多溴二苯醚 PBDE	
包装 Packaging	O	O	O	O	O	O	
塑料外壳 / 组件 Plastic housing / parts	O	O	O	O	O	O	
真空油 Vacuum oil	O	O	O	O	O	O	
电池 Battery	O	O	O	O	O	O	
玻璃 Glass	O	O	O	O	O	O	
电子电气组件 Electrical and electronic parts	X	X	X	O	O	O	
控制器 / 测量设备 Controller / measuring device	X	O	X	O	O	O	
金属外壳 / 组件 Metal housing / parts	X	O	O	O	O	O	
电机 Motor	X	O	O	O	O	O	
配件 Accessories	X	O	O	O	O	O	

The logo for vacuubrand, featuring the word "vacuubrand" in a lowercase, sans-serif font, enclosed within a blue oval border.

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

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O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.

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X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by VACUUBRAND may enter into further devices (e.g., rotary evaporator) or can be used together with other appliances (e.g., usage as booster pumps).

With these products and appliances in particular, please note the EFUP labeled on these products. VACUUBRAND will not take responsibility for the EFUP of those products and appliances.

Place. date: Wertheim, 04/24/2017

A handwritten signature in black ink, appearing to read "F. Gitmans".

(Dr. F. Gitmans)
Managing Director

ppa.

A handwritten signature in black ink, appearing to read "J. Dirscherl".

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Disclaimer: Our technical literature is only intended to inform our customer. The validity of general empirical values and results obtained under test conditions for specific applications depend on a number of factors beyond our control. It is therefore strictly the users' responsibility to very carefully check the validity of application to their specific requirements. No claims arising from the information provided in this literature will, consequently, be entertained.

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VACUUBRAND GMBH + CO KG
-Technology for Vacuum Systems-

© 2018 VACUUBRAND GMBH + CO KG Printed in Germany

Alfred-Zippe-Str. 4 · 97877 Wertheim / Germany
T +49 9342 808-0 · F +49 9342 808-5555
info@vacuubrand.com · www.vacuubrand.com

20901409 / 01/07/2018

Documents are only to be used and distributed completely and unchanged. It is strictly the users' responsibility to check carefully the validity of this document with respect to his product.

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VACUU•LAN® (US-Reg.No 3,704,401), VACUU•BUS®, VACUU•CONTROL®, Peltronic®, VARIO® (US-Reg.No 3,833,788), VACUUBRAND® (US-Reg.No 3,733,388) and also the shown company logos are registered trademarks of VACUUBRAND GMBH + CO KG in Germany and/or other countries.