# VAKUUM-SCHLAUCH VACUUM HOSE – TUYAU DE VIDE

*PTFE-AS*

**Hinweise zu Installation und Betrieb**  
**Notes concerning installation and use**  
**Notes concernant l'installation et l'utilisation**

# Vakuum-Schlauch PTFE-AS

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# Vacuum hose PTFE-AS

## Notes concerning installation and use

PTFE-AS is antistatic according to BS 5958:1991 / EN ISO 8031.

When installing the vacuum hose do not under-run the minimum bending radius of the hose (see "Technical data"). Do not kink. Design the hose length appropriately and bear in mind the own weight of the hose.

Avoid traction and abrasion by varying the hose length or the mounting, by using special armatures or by supporting the vacuum hose appropriately.

Install hose without twisting, no torsional stress must occur.

After installation check hose regularly for mechanical damages. Do not use damaged hoses, replace with undamaged ones.

If the hose is to be used in potentially explosive atmospheres, ensure that an antistatic conductive connection to parts connected to ground is established. The vacuum hose must not be used to connect other devices to ground. If necessary, those devices have to be connected to ground separately.

## Technical data

Vacuum hose PTFE-AS		NW 16	NW 25
Nominal inner diameter	mm	16	25
Length	mm	500 or 1000 or 2500	
Minimum bending radius dynamic <sup>1,2</sup>	mm	300	400
Minimum bending radius static <sup>1</sup>	mm	150	200
Connections		small flange DN 16	small flange DN 25
Permitted range of ambient temperature	°C	-10 to +40	
Permitted temperature range of pumped media	°C	- 50 to +100	
Leak rate	mbar l/s	< 10 <sup>-3</sup>	
Permitted pressure range during operation		vacuum to max. 1.2 bar absolute	

1: Avoid torsional stress (twisting) of the hose. Do not kink, damage or mechanically squeeze the vacuum hose.

2: Hose must not chafe (abrasion).

Component	Wetted material
Hose	PTFE antistatic
Flange connection	stainless steel
Retaining ring	stainless steel (not wetted)

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