ELEKTROMAGNETVENTIL SOLENOID OPERATED VALVE ÉLECTROVANNE

VCL - C3

Instructions for use

Inhalt / Contents / Contenu

Unbedingt beachten!	3
Technische Daten	4
Safety information!	5
Technical data	6
Notes importantes!	7
Données techniques	8
China RoHS	9

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Safety information!

NOTICE

Use the solenoid operated in-line valve for vacuum control by connecting it to a vacuum controller.

- ☞ Connect the solenoid operated valve only to a suitable vacuum controller (see "Technical data").

The solenoid operated valve is designed to be connected to the VACUUBRAND BUS system VACUU·BUS.

Remove all packing material from the packing box. Remove the product from its packing-box and retain all packaging for inspection until the equipment is inspected and tested. Inspect the equipment promptly and carefully. If the equipment is damaged, notify the supplier and the carrier in writing within three days. State the item number of the product together with the order number and the supplier's invoice number.

Use the equipment **for the intended use only**, i.e. as in-line valve in VACUU·LAN® Local Area Vacuum Networks or in vacuum systems. Do not use at overpressure. Do not use to distribute fluids or solids!

Comply with all applicable and relevant safety requirements (regulations and guidelines). **Implement the required actions and adopt suitable safety measures.**

Ensure that **maintenance** is done only by suitably trained and supervised technicians. As laid down in the statutory regulations (occupational, health and safety regulations and regulations for environmental protection), equipment returned to the manufacturer can be repaired only under certain conditions.

! DANGER

- ➔ Adopt suitable measures to prevent the release of **dangerous, toxic, explosive, corrosive, noxious or polluting fluids, vapours and gases**.
- ➔ **Attention:** The equipment may be contaminated with chemicals, which have been processed during operation. Ensure that the equipment is completely decontaminated before maintenance commences. Take adequate precautions to protect people from the effects of dangerous substances if contamination has occurred. Use appropriate protective clothing, safety goggles and protective gloves.

! WARNING

- ☞ Use only **OEM spare parts and accessories**. Otherwise safety and performance of the equipment might be reduced.
- ☞ **Do not use the equipment if it is damaged.**
- ☞ Check that the system is mechanically stable and that all fittings are secure.
- ☞ Do not disassemble any valve if under vacuum!

! CAUTION

- Check that the power source and the equipment match in voltage, phase, and frequency (see "Technical data").
- Ensure that the materials of the wetted parts are compatible with the media, see "Technical data".
- **If idle, the solenoid operated valve is closed.**
The valve is open if energized and if Bus signal is applied.
- Comply with max. permissible temperatures and pressures (see "Technical data").
- Check regularly that the connection between valve and vacuum system is leak-tight, especially prior to initial operation.

Assembly:

- ➔ Do not insert a sealing ring or a non-return valve.
- ➔ Tighten the **union nut** hand-tight; retighten after some days if necessary.

Technical data

Type	20636668
Maximum permissible range of supply voltage	24 VDC $\pm 10\%$
Power draw (operation)	6 W
Current draw	approx. 0.22 A
Electrical connection	connector for VACUUBRAND BUS system VACUU·BUS
Operating cycles per minute	maximum 50
Valve currentless	valve closed
Nominal diameter	DN 6
Max. permissible operating pressure	1.5 bar
Connection	union nut M35 x 1.5
Max. permissible ambient temperature	40 °C
Max. permissible gas temperature	80 °C
Degree of protection IEC 529	IP 65
Leak rate	$1 \cdot 10^{-2}$ mbar·l/s
Dimensions L x W x H (without cable)	60 mm x 41 mm x 68 mm
Cable length	2 m
Weight (without cable)	0.21 kg

Components	Wetted materials
Valve body	PVDF
Seal material	PTFE / fluorine elastomer

We reserve the right for technical modifications without prior notice!

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DECLARATION OF CONFORMITY – China RoHS 2

VACUUBRAND GMBH + CO KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a “Product Conformity Assessment” (PCA) procedure was performed. As defined in GB 26572-2025 the “Maximum Concentration Value” limits (MCV) apply to these restricted substances:

• Lead (Pb):	0.1%
• Mercury (Hg):	0.1%
• Cadmium (Cd):	0.01%
• Hexavalent chromium (Cr(+VI)):	0.1%
• Polybrominated biphenyls (PBB):	0.1%
• Polybrominated diphenyl ether (PBDE):	0.1%
• Dibutyl phthalate (DBP):	0.1%
• Benzyl butyl phthalate (BBP):	0.1%
• Bis(2-ethylhexyl) phthalate (DEHP):	0.1%
• Diisobutyl phthalate (DIBP):	0.1%

Environmentally Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmentally Friendly Use Period for VACUUBRAND products is 40 years.

VACUUBRAND®

部件名称 Part name	有害物质 / Hazardous substances									
	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr ⁶⁺)	多溴联苯 (PBB)	多溴二苯醚 (PBDE)	邻苯二甲酸二(2-乙基己)酯 (DEHP)	邻苯二甲酸丁苄酯 (BBP)	邻苯二甲酸二丁酯 (DBP)	邻苯二甲酸二异丁酯 (DIBP)
包装 Packaging	O	O	O	O	O	O	O	O	O	O
塑料外壳 / 组件 Plastic housing / parts	O	O	O	O	O	O	O	O	O	O
真空油 Vacuum oil	O	O	O	O	O	O	O	O	O	O
电池 Battery	O	O	O	O	O	O	O	O	O	O
玻璃 Glass	X	O	O	O	O	O	O	O	O	O
电子电气组件 Electrical and electronic parts	X	O	O	O	O	O	O	O	O	O
控制器 / 测量设备 Controller / measuring device	X	O	O	O	O	O	O	O	O	O
金属外壳 / 组件 Metal housing / parts	X	O	O	O	O	O	O	O	O	O
电机 Motor	X	O	O	O	O	O	O	O	O	O
配件 Accessories	X	O	O	O	O	O	O	O	O	O

- O: 表示该有害物质在该部件所有均质材料中的含量均在 GB 26572-2025 标准规定的限量要求以下。
O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB 26572-2025.
- X: 表示该有害物质在该部件至少一个均质材料中的含量超出 GB 26572-2025 标准规定的限量要求。
X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB 26572-2025.

Declaration of Conformity – China RoHS 2**V5_November 2025**

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Batteries, glassware and accessories might not be content of the enclosed device and may have its own EFUP-marking and/or might be maintaining parts with changing EFUP-marking.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), polybrominated diphenyl ethers (PBDE), dibutyl phthalate (DBP), benzyl butyl phthalate (BBP), bis(2-ethylhexyl) phthalate (DEHP), Diisobutyl phthalate (DIBP).

Products manufactured by VACUUBRAND may enter into further devices (e.g., rotary evaporator) or can be used together with other appliances (e.g., usage as booster pumps).

With these products and appliances in particular, please note the EFUP labeled on these products.


VACUUBRAND will not take responsibility for the EFUP of those products and appliances.

Place, date: Wertheim, 14.11.2025

(Dr. Constantin Schöler)
Managing Director

ppa.


(Jens Kaibel)
Technical Director

VACUUBRAND GMBH + CO KG

Alfred-Zippe-Str. 4
97877 Wertheim
Germany

Tel.: +49 9342 808-0
Fax: +49 9342 808-5555
E-Mail: info@vacuubrand.com
Web: www.vacuubrand.com

Manufacturer:

VACUUBRAND GMBH + CO KG
Alfred-Zippe-Str. 4
97877 Wertheim
GERMANY

Phone:

Head office +49 9342 808-0
Sales +49 9342 808-5550
Service +49 9342 808-5660

Fax: +49 9342 808-5555

Email: info@vacuubrand.com

Web: www.vacuubrand.com