

VacuHandControl VHC^{pro}

[Lien vers le site web](#)



Description du produit

La poignée d'aspiration VHC^{pro} est légère, ergonomique et autoclavable à 121°C. Elle est livrée prête à l'emploi avec 2,5 m de tuyau silicone très flexible, des adaptateurs pour cônes, pipettes graduées et pipettes Pasteur, un adaptateur mural pour reposer la poignée (pour maintenir le poste de travail libre) et des clips adhésifs pour guider le tuyau d'aspiration dans les postes de sécurité microbiologiques. Le mécanisme de la poignée VHC n'est pas en contact avec les fluides contaminés, il est fiable et exempt de fuites.

Caractéristiques

- léger et ergonomique pour une manipulation confortable
- position d'aspiration permanente pour un débit continu
- aspiration via tube traversant pour manipulation sans fuite
- autoclavable à 121°C et 2 bar

Référence de commande

Désignation	Tension fréquence nominales	Type de fiche secteur	Référence
VHC ^{PRO}			20688061

Livraison

Poignée d'aspiration ergonomique avec adaptateur pour différentes pipettes et pointes, levier pour aspiration manuelle et bouton pour aspiration en continu, tubulure de 2,5 m.

Accessoires

Caractéristiques techniques

Température ambiante maxi admissible (service)	10 - 40 °C
Température ambiante maxi admissible (stockage)	-10 - 60 °C
Pression max. à l'échappement (abs.)	1.1 bar
Raccord côté aspiration	Poignée d'aspiration
Raccord côté pression	--
Classe de protection IP selon IEC 60529	IP 20
Certification NRTL	Non

Plan de masse



VACUUBRAND GMBH + CO KG
P.O. Box 1664
97877 Wertheim
T +49 9342 808 5550
F +49 9342 808 5555

info@vacuubrand.com
www.vacuubrand.com

VACUUBRAND®

VACUUBRAND is part of the Brand group, a manufacturer of high-quality and innovative laboratory equipment and vacuum technology with a global team of approx. 1,000 employees.

With highly motivated and qualified employees, we manufacture the world's most comprehensive product family for vacuum generation, measurement and control for the rough and fine vacuum range in the laboratory.

All with one goal:

Enabling our customers to run reliable, safe and efficient processes. This is why we commit ourselves to continuous innovation, first-class performance, and outstanding quality.

Engineering excellence is central to our ability to deliver products and solutions renowned for their quality, durability, and reliability. Innovative technology, high quality materials and extensive testing deliver unmatched performance - engineered and manufactured in Germany.

Together with you, our experts will find the best solution for your application and provide you with first-class service along the way.

BRAND (Shanghai) Trading Co., Ltd.
Shanghai, China

T +86 21 6422 2318
info@brand.com.cn
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.
Mumbai, India

T +91 22 42957790
customersupport@brand.co.in
www.brand.co.in

BRANDTECH® Scientific, Inc.
Essex, CT. United States of America

T +1 860 767 2562
info@brandtech.com
www.brandtech.com