

Capteur de vide VACUU·SELECT Sensor

[Lien vers le site web](#)

Description du produit

Référence de commande

Désignation	Tension fréquence nominales	Type de fiche secteur	Référence
VACUU·SELECT Sensor	VACUU·BUS DC	VACUU·BUS®	20700020

Livraison

Vacuum gauge head VACUU·SELECT Sensor with venting valve.

Caractéristiques techniques

Plage de mesure (mbar)	1080 - 0.1 mbar
Plage de mesure (torr)	810 - 0.1 torr
Principe de mesure	Membrane en céramique d'alumine, capacitif, ind. du gaz, mesure absolue
Précision de mesure	< +- 1 mbar/hPa/torr / +- 1 chiffre (après ajustage, à temp. constante)
Dérive de température	< 0.15 mbar/hPa/0.11 torr /K
Raccord de commande	Petide bride KF DN 16, PTFE raccord tuyau 10/8 mm, embout DN 6/10 mm
Température ambiante maxi admissible (service)	10 - 40 °C
Température ambiante maxi admissible (stockage)	-10 - 60 °C
température du milieu gazeux (fonctionnement en continu)	40 °C
Temp. Max. fluide pompé temps court	80 °C
Matériau du boîtier	Boîtier en plastique stable avec bonne résistance aux produits chimiques
Classe de protection IP selon IEC 60529	IP 41
Classe de protection selon UL50E	Type 2
Poids	0.15 kg
Longueur du câble de connexion	2 m
Conformité ATEX	II 3/- G Ex h IIC T4 Gc X Internal Atm. only
Certification NRTL	Oui
Interface	VACUU-BUS

Plan de masse



VACUUBRAND GMBH + CO KG
P.O. Box 1664
97877 Wertheim
T +49 9342 808 5550
F +49 9342 808 5555

info@vacuubrand.com
www.vacuubrand.com

VACUUBRAND®

VACUUBRAND is part of the Brand group, a manufacturer of high-quality and innovative laboratory equipment and vacuum technology with a global team of approx. 1,000 employees.

With highly motivated and qualified employees, we manufacture the world's most comprehensive product family for vacuum generation, measurement and control for the rough and fine vacuum range in the laboratory.

All with one goal:

Enabling our customers to run reliable, safe and efficient processes. This is why we commit ourselves to continuous innovation, first-class performance, and outstanding quality.

Engineering excellence is central to our ability to deliver products and solutions renowned for their quality, durability, and reliability. Innovative technology, high quality materials and extensive testing deliver unmatched performance - engineered and manufactured in Germany.

Together with you, our experts will find the best solution for your application and provide you with first-class service along the way.

© VACUUBRAND GMBH + CO KG

BRAND (Shanghai) Trading Co., Ltd.
Shanghai, China

T +86 21 6422 2318
info@brand.com.cn
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.
Mumbai, India

T +91 22 42957790
customersupport@brand.co.in
www.brand.co.in

BRANDTECH® Scientific, Inc.
Essex, CT. United States of America

T +1 860 767 2562
info@brandtech.com
www.brandtech.com