

# VacuHandControl

## VHC<sup>pro</sup>

[网站链接](#)

### 产品说明

VHC<sup>pro</sup> 是一款按人体工程学设计的轻质吸液手柄，可进行121 °C高压灭菌。整套手柄包含2.5米长流经液体的柔软的硅胶管，可连用刻度移液管、吸管和巴斯德管的适配器，吸头适配器，壁挂支架（无需占据台面）及带粘性背胶的夹子（可将吸液软管固定以保持生物安全柜的整洁）。VHC<sup>pro</sup>的机理使手柄柄身不接触废液，使用更灵活，完全不用担心泄露问题。

### 性能特点

- 重量轻、人体工程学设计，操作非常便捷
- 持续吸液功能，可大量吸液
- 废液流经管路，完全不用担心漏液
- 高压灭菌，121 °C和2 bars

### 订购信息

产品名称	额定电源电压 工作频率	电源插头类型	订货号
VHC <sup>pro</sup>			20688061

### 货品配置

人体工程学设计手柄，带吸头接头，可以适配一系列普通吸头，可进行手动吸液控制，或旋转钮至持续吸液功能，带2.5米软管

### 配件

## 技术指标

的允许环境温度范围 (运行)	10 - 40 °C
的允许环境温度范围 (存放)	-10 - 60 °C
最大出气口压力 (绝对压力)	1.1 bar
进气口连接	手柄
出气口连接	--
防护等级	IP 20
NRTL 标准	否

尺寸规格表



VACUUBRAND GMBH + CO KG  
P.O. Box 1664  
97877 Wertheim  
T +49 9342 808 5550  
F +49 9342 808 5555

info@vacuubrand.com  
www.vacuubrand.com

**VACUUBRAND®**

VACUUBRAND is part of the Brand group, a manufacturer of high-quality and innovative laboratory equipment and vacuum technology with a global team of approx. 1,000 employees.

With highly motivated and qualified employees, we manufacture the world's most comprehensive product family for vacuum generation, measurement and control for the rough and fine vacuum range in the laboratory.

**All with one goal:**

Enabling our customers to run reliable, safe and efficient processes. This is why we commit ourselves to continuous innovation, first-class performance, and outstanding quality.

Engineering excellence is central to our ability to deliver products and solutions renowned for their quality, durability, and reliability. Innovative technology, high quality materials and extensive testing deliver unmatched performance - engineered and manufactured in Germany.

Together with you, our experts will find the best solution for your application and provide you with first-class service along the way.

BRAND (Shanghai) Trading Co., Ltd.  
Shanghai, China

T +86 21 6422 2318  
info@brand.com.cn  
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.  
Mumbai, India

T +91 22 42957790  
customersupport@brand.co.in  
www.brand.co.in

BRANDTECH® Scientific, Inc.  
Essex, CT. United States of America

T +1 860 767 2562  
info@brandtech.com  
www.brandtech.com